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

Analitik Evre


FIGURE 6-2 Process flow for the three phases of blood gas and analyte measurement.

pH arterial	Tam kan	Genel	≤ 7.200	≥ 7.600	pH
pCO ₂ , arterial	Tam kan	Genel	≤ 20.0	≥ 70.0	mmHg
pO ₂ , arterial	Tam kan	Genel	≤ 40.0	-	mmHg

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Delost, M. (2014). Blood gas and critical care analyte analysis. *Equipment for Respiratory Care; Volsko, TA, Chatburn, RL, El-Khatib, MF, Eds.*
<https://dosyamerkez.saglik.gov.tr/Eklenti/15356/0/karar-siniri-esik-deger-kritik-deger-panik-deger-ve-olcum-birimlerinin-harmonizasyonu-proseduru.pdf> (Accessed 1 November 2024).

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Mathew, B., & Lakshminrusimha, S. (2016). Noninvasive monitoring of gas exchange. *Assisted Ventilation of the Neonate. An Evidence-Based Approach to Newborn Respiratory Care*, JP Goldsmith, EH Karotkin, M. Keszler, and GK Suresh, Eds, 97-107.

Davis, M. D., Walsh, B. K., Sittig, S. E., & Restrepo, R. D. (2013). AARC clinical practice guideline: blood gas analysis and hemoximetry: 2013. *Respiratory care*, 58(10), 1694-1703.

Table 1. Test parameters and methods on blood gas analyzers.

Description	Parameter	Method
Traditional blood gas		
	pO ₂	Clark Electrode
	pCO ₂	Severinghaus electrode
	pH	Sanz electrode
Electrolytes		
	Na ⁺	Ion selective electrode
	K ⁺	Ion selective electrode
	Cl ⁻	Ion selective electrode
	iCa ⁺⁺	Ion selective electrode
Metabolites		
	Glucose	Immobilized enzyme/electrode
	Lactate	Immobilized enzyme/electrode
	Urea nitrogen	Immobilized enzyme/electrode
	Creatinine	Immobilized enzyme/electrode
Co-oximetry		
	Total hemoglobin	Differential spectrophotometry
	FO ₂ Hb	Differential spectrophotometry
	COHB ^a	Differential spectrophotometry
	MetHb ^b	Differential spectrophotometry
Calculated parameters		
	HCO ₃ ⁻	From pH and pCO ₂
	Anion gap	From Na, Cl, and calculated HCO ₃ ⁻
	"O ₂ Sat"	Estimated from pO ₂ and normal oxyhemoglobin dissociation curve

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Severinghaus, J. W. (2002). The invention and development of blood gas analysis apparatus. *The Journal of the American Society of Anesthesiologists*, 97(1), 253-256.

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Table C1. Measurement Technologies Used in Instruments for Blood Gases, pH, and Related Analytes

	Electrochemical			Optical
	Potentiometry	Amperometry	Conductometry	Absorption Spectrometry
pH and Gases	pH, $p\text{CO}_2$	$p\text{O}_2$		
Electrolytes	Na^+ , K^+ , Cl^- , Ca^{++} , Mg^{++}			
Metabolites		glucose, lactate		
Hematocrit, Hemoglobin, Hemoglobin Fractions			hematocrit	hemoglobin hemoglobin fractions

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
Siemens Healthineers Inc., Germany

Company	Product	Testing analytes	Measured parameters	Calculated parameters	Sample type	Sample size
Siemens Healthineers Inc., Germany	RAPIDLab 1240	Blood gas	—	—	Heparinized whole blood	90–175 µL
	RAPIDLab1245	Blood gas CO-oximetry Neonatal total bilirubin	—	—	Heparinized whole blood	
	RAPIDLab 1260	Blood gas Electrolytes Metabolites	—	—	Heparinized whole blood	
	RAPIDLab 1265	Blood gas Electrolytes Metabolites CO-oximetry Neonatal total bilirubin	—	—	Heparinized whole blood	
	RAPIDLab 348EX System	Blood gas pH CO-oximetry Metabolites Hematocrit	—	—	Whole blood Dialysate fluid	50–95 µL
	epoc Blood Analysis System (Formerly under Alere)	Blood gas Electrolytes Metabolites	pH, pCO ₂ , pO ₂ , TCO ₂ , Na ⁺ , K ⁺ , Ca ⁺⁺ , Hct, Glu, Lac, Crea, Urea	Agap, AGapK, CHCO ₃ , cTCO ₂ , BE(ecf), BE(b), cSO ₂ , cHgb, eGFR, eGFR-a	Whole blood	92 µL

Dedeoglu, A., Kaya, S. I., Bakirhan, N. K., & Ozkan, S. A. (2020). Nanotechnological approaches and materials in commercial biosensors. In *Commercial biosensors and their applications* (pp. 301-353). Elsevier. RAPIDLab 1200 Systems - Siemens Healthineers <https://www.siemens-healthineers.com/tr/blood-gas/blood-gas-systems/rapidlab-1200-systems> (Accessed 1 November 2024).

Instrumentation Laboratory Co.,USA

Company	Product	Testing analytes	Measured parameters	Calculated parameters	Sample type	Sample size
Instrumentation Laboratory Co., USA	GEM Premier 3000	Blood gas pH CO-oximetry Metabolites Hematocrit	pH, pCO ₂ , pO ₂ , Na ⁺ , K ⁺ , Ca ⁺⁺ , Glu, Lac, Hct	HCO ₃ ⁻ , HCO ₃ ⁻ std, TCO ₂ , BE(B) in vitro, BE(ecf) in vivo, SO ₂ c, Ca ⁺⁺ (7.4), THbc, A-aDO ₂ , pAO ₂ , paO ₂ /pAO ₂ , RI, CaO ₂ , CvO ₂ , CcO ₂ , a-vDO ₂ , Qsp/Qt, P50	Whole blood	—
	GEM Premier 4000	pH Blood gas Electrolytes Metabolites CO-oximetry	pH, pCO ₂ , pO ₂ , Hct, Na, K, Cl, iCa, lactate, glucose, tHb(c), O ₂ Hb, COHb, MetHb, HHb, tBili	TCO ₂ , BEecf (in vivo), BE(B) (in vivo), tHb(c), Ca ⁺⁺ (7.4), anion gap, P/F ratio, pAO ₂ , CaO ₂ , CvO ₂ , P50, O ₂ ct, O ₂ cap, SO ₂ , SO ₂ (c), HCO ₃ -std, HCO ₃ -(c), A-aDO ₂ , paO ₂ /pAO ₂ , RI, CcO ₂ , a-vDO ₂ , Qsp/Qt(est), Qsp/Qt, Hct(c)	Whole blood	65 µL
	GEM Premier 5000	pH Blood gas Electrolytes Metabolites CO-oximetry	pH, pCO ₂ , pO ₂ , Na ⁺ , K ⁺ , Ca ⁺⁺ , Cl ⁻ , Glu, Lac, Hct, tHb, O ₂ Hb, COHb, MetHb, HHb, tBili, SO ₂	BE(B), BE(ecf), tHb(c), Ca ⁺⁺ (7.4), anion gap (AG), P/F ratio, pAO ₂ , CaO ₂ , CvO ₂ , p50, O ₂ cap, SO ₂ (c), O ₂ ct, HCO ₃ - std, TCO ₂ , HCO ₃ -(c), A-a DO ₂ , paO ₂ /pAO ₂ , RI, CcO ₂ , a-vDO ₂ , Qsp/Qt (est), Qsp/Qt, Hct(c)	Whole blood	
	GEM Premier 3500	pH Blood gas Electrolytes Metabolites	pH, pCO ₂ , pO ₂ , Na ⁺ , K ⁺ , Ca ⁺⁺ , Glu, Lac, Hct	HCO ₃ ⁻ , HCO ₃ ⁻ std, TCO ₂ , BE(B) in vitro, BE(ecf) in vivo, SO ₂ c, Ca ⁺⁺	Whole blood	

GEM Premier 3500 j Instrumentation Laboratory Worldwide, <https://www.werfen.com/na/en/point-of-care-testing-devices/gem-premier-3500-poc-blood-gas-analyzer> (Accessed 1 November 2024).

Dedeoglu, A., Kaya, S. I., Bakirhan, N. K., & Ozkan, S. A. (2020). Nanotechnological approaches and materials in commercial biosensors. In *Commercial biosensors and their applications* (pp. 301-353). Elsevier.

Three-year monitoring and comparison of results from two identical blood gas analyzers

Yun Huang^{a,b,*}, Robert Dean^a, Yvonne Dubbelman^a, Anne Vincent^a,
Ying-pui Michael Chan^b

^a Clinical Laboratories, Kingston Health Sciences Centre, 76 Stuart Street, Kingston, ON, Canada

^b Department of Pathology and Molecular Medicine, Queen's University, 76 Stuart Street, Kingston, ON, Canada

ARTICLE INFO

Keywords:

Blood gas analyzer
Result interchangeability
Between-analyzer comparison
Point-of-care testing

ABSTRACT

Objectives: Measurement comparability between blood gas analyzers within a laboratory is of utmost importance. This study analyzed the data obtained from a three-year period.

Design and methods: For quality monitoring one blood sample was tested on two identical blood gas analyzers at each of three shifts/day for three years. Deming regression analysis was used to determine result correlation and statistical identity between the two analyzers for pH, pCO₂, pO₂, sodium, potassium, chloride, ionized calcium, glucose, and lactate. Failures in the two-analyzer comparison were determined by the performance limits from the Institute of Quality Management in Healthcare (IQMH) and from the manufacturer respectively.

Results: Correlation coefficients were greater than 0.96 (0.9622–0.9975) for all tested analytes. The measurements of every analyte on both analyzers were statistically identical. In the two-analyzer comparison failure numbers/1000 tests for pO₂ and glucose varied with the performance limits (IQMH: 0.6 and 49.2; the manufacturer: 19.3 and 4.4, respectively). In addition, persistent glucose failures (>5/week) between the two analyzers occurred occasionally.

Conclusions: Results of all tested analytes between the two blood gas analyzers were interchangeable. Recurring glucose discrepancies might be a result of different lots of cartridges used on each analyzer, which were not identified during the initial installation.

Radiometer Inc., USA

Company	Product	Testing analytes	Measured parameters	Calculated parameters	Sample type	Sample size
	ABL800 FLEX Analyser	—	pH, pCO ₂ , pO ₂ , SO ₂ , ctHb, FO ₂ Hb, FCOHb, FMetHb, FHHb, FHbF, cK ⁺ , cNa ⁺ , cCa ²⁺ , cCl ⁻ , cGlu, cLac, cCrea, ctBil	Anion gap, anion Gap(K ⁻), cCa ²⁺ (pH = 7.40), pO ₂ (a,T)/FO ₂ (I), pO ₂ (a)/FO ₂ (I), pO ₂ (a/A,T), pO ₂ (a/A), pO ₂ (A-a), pO ₂ (A-a), p50(st)	Whole blood	35 µL
	ABL90 FLEX PLUS Analyser	—	pH, pCO ₂ , pO ₂ , SO ₂ , ctHb, FO ₂ Hb, FCOHb, FMetHb, FHHb, FHbF, cK ⁺ , cNa ⁺ , cCa ²⁺ , cCl ⁻ , cGlu, cLac, cCrea, cUrea, ctBil	—	—	65 µL
Radiometer Inc., USA	ABL9 Blood Gas Analyser	—	pH, pCO ₂ , pO ₂ , Na ⁺ , K ⁺ , Cl ⁻ , iCa, Glu, BUN, Hct, creatinine, lactate	CH ⁺ , ctHb, CHCO ₃ (P), CHCO ₃ (P,st), cBase(B), cBase(Ect), cBase(B,ox), cBase(F,ox), cCa ²⁺ (7.40), etCO ₂ (B), etCO ₂ (P), anion gap, anion gap (K ⁻), Measured parameters, pO ₂ (A), pO ₂ (A-a), pO ₂ (a/A), pO ₂ (a)/FO ₂ I, SO ₂ , etO ₂ , RI	Whole blood	70 µL

Dedeoglu, A., Kaya, S. I., Bakirhan, N. K., & Ozkan, S. A. (2020). Nanotechnological approaches and materials in commercial biosensors. In *Commercial biosensors and their applications* (pp. 301-353). Elsevier. ABL800 FLEX blood gas analyzer <https://www.radiometer.com/en/products/blood-gas-testing/abl800-flex-blood-gas-analyzer> (Accessed 1 November 2024).

Abbott i-STAT

Nova Biomedical, Waltham, MA

4 Equipment

- 1 i-STAT ALINITY INSTRUMENT: Used to perform cartridge testing, reviewing test results, and conducting quality control (QC) testing.
 - 2 i-STAT ALINITY BASE STATION: Used to recharge the battery installed in the i-STAT Alinity.
 - 3 i-STAT CARTRIDGES: Contains sensors and reagents for all patient and quality testing.
 - 4 i-STAT ALINITY RECHARGEABLE BATTERY: Provides main power source to the instrument.
 - 5 i-STAT ALINITY ELECTRONIC SIMULATOR: Provides an independent check on the instrument's thermal controls and success of software updates.
 - 6 i-STAT ALINITY PORTABLE PRINTER: Used to print records stored in the instrument.
- G3+ cartridges: 3P7825 blood gas for respiratory team
 - EG7+ cartridges: 3P7625 blood gas and electrolyte
 - CG4+ cartridges: 3P8525 blood gas and lactate
 - iSTAT printer paper (6F17-11) is obtained from POCT
 - Abbott Tricontrol Level 1 (5P7101) and Level 3 (5P7301)

SYSTEM COMPONENTS

Füzéry AK, Elian FA, Kost GJ. A review of temperature-related challenges and solutions for the Abbott i-STAT and Siemens Healthineers epoc devices. Clin Biochem. 2023 May;115:49-66. doi: 10.1016/j.clinbiochem.2022.08.013. Epub 2022 Sep 5. PMID: 36067872.

<https://www.globalpointofcare.abbott/us/en/product-details/apoc/istat-eg7plus-us.html>
(Accessed 1 November 2024).

<https://www.novabiomedical.com/press/100323.php>
(Accessed 1 November 2024).


Roche, Roche Diagnostics GmbH, Mannheim, Germany

- **Cobas® b 221 system**
- **Cobas® b 123 POC system**
- pH
- Blood Gas (BG): PO₂, PCO₂
- Electrolyte (ISE): Na⁺, K⁺, Cl⁻, iCa²⁺
- Hematocrit (Hct)
- Metabolites (MSS): Glucose, Lactate, Urea/BUN (only on **cobas b 221<6>** system)
- Total Hemoglobin (tHb)
- Oxygen saturation (SO₂)
- Hemoglobin derivative COOX (O₂Hb, HHb, COHb, MetHb)
- Bilirubin (neonatal)

Roche Diagnostics, 2022. cobas b 221 system operators manual version 19.0. pp. 35.

<https://diagnostics.roche.com/global/en/products/instruments/cobas-b-221-2-system-ins-391.html> (Accessed 1 November 2024).

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Srinivasa N. Raja, M.D., Editor

Anesthesiology 2002; 97:253-6

© 2002 American Society of Anesthesiologists, Inc. Lippincott Williams & Wilkins, Inc.

The Invention and Development of Blood Gas Analysis Apparatus

John W. Severinghaus, M.D., F.R.C.A., Dr.Med.H.C.*

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Fig. 1. Leland Clark with his invention of a method of defoaming in a bubble oxygenator (1952).

Fig. 2. First three-function blood gas apparatus (pH, partial pressure of carbon dioxide [Pco₂], and partial pressure of oxygen [Po₂] electrodes) built by Severinghaus and Bradley (1959).

References

1. Osborn JJ: Experimental hypothermia: Respiratory and blood pH changes in relation to cardiac function. *Am J Physiol* 1953; 175:389-98
2. Severinghaus JW, Stupfel M, Bradley AF Jr: Accuracy of pH and Pco₂ determinations. *J Appl Physiol* 1956; 9:189-96
3. Stow RW, Randall BF: Electrical measurement of the Pco₂ of blood (abstract). *Am J Physiol* 1954; 179:678
4. Stow RW, Baer RF, Randall B: Rapid measurement of the tension of carbon dioxide in blood. *Arch Phys Med Rehabil* 1957; 38:646-50
5. Clark LC Jr: Monitor and control of blood and tissue O₂ tensions. *Trans Am Soc Artif Intern Organs* 1956; 2:41-8
6. Severinghaus JW, Bradley AF Jr: Electrodes for blood pO₂ and pCO₂ determination. *J Appl Physiol* 1958; 13:515-20
7. Astrup P: A simple electrometric technique for the determination of carbon dioxide tension in blood and plasma, total content of carbon dioxide in plasma and bicarbonate content in "separated" plasma at a fixed carbon dioxide tension. *Scand J Clin Lab Invest* 1956; 8:33-43
8. Siggaard-Andersen O, Engel K, Jorgensen K, Astrup P: A micro method for determination of pH, carbon dioxide tension, base excess and standard bicarbonate in capillary blood. *Scand J Clin Lab Invest* 1960; 12:172-6
9. Severinghaus JW: Blood gas calculator. *J Appl Physiol* 1966; 21:1108-16
10. Roughton FJW, Severinghaus JW: Accurate determination of O₂ dissociation curve of human blood above 98.7% saturation with data on O₂ solubility in unmodified human blood from 0° to 37°C. *J Appl Physiol* 1973; 35:861-9
11. Severinghaus JW: Simple, accurate equations for human blood O₂ dissociation computations. *J Appl Physiol* 1979; 46:599-602